

BBS Revolution™

The Next Generation Bladder Scanner

Fully automated bladder scanner with a unique, user friendly design.


- Automatically detects male/female anatomy
- Automatically locates bladder without aiming
- Quickly provides results you can trust

Improved Technology + Reduced Total Cost of Ownership +
Decreased Infection Risk
= BETTER SOLUTION



BBS Revolution Clinical Benefits

- ...Saves nurses time and frustration from setup through cleanup
- ...Reduces patient discomfort due to probe design and quick test result
- ...Simplifies disinfection requirements via the wireless operation
- ...Increases compliance with infection control protocols
- ...Improves accuracy and uptime with self-calibration mode
- ...Minimizes user training through its elegantly simple design



Waging War Against Healthcare-Associated Infections

- ✓ Lower Costs
- ✓ Better Outcomes
- ✓ More satisfied patients and providers

dBMED_x

BBS Revolution™

The Next Generation Bladder Scanner

Wireless
Scanner

8 transducers,
180 degree view angle

Internal Battery,
lasts up to 48 hours of use
or 500 scans




It takes just one quick scan to get an accurate, reliable volume measurement.

Follow these simple steps:

1. Apply ample gel on the abdomen
2. Place the probe firmly against the lower abdomen
3. Press the button to take a reading

IT IS THAT EASY!



There is no need to select modes or aim the probe, it is ONE AND DONE™ with the BBS Revolution™.

CONTACT US FOR A DEMO

- www.dBMEDx.com
- 720-515-6630
- info@dbmedx.com

